

CERTIFICATE OF ANALYSIS

Product Human Serum Lipid Stripped - Mixed Gender, Pooled
Item No. 175CM-100MLX
Lot No. 175CM-XXA24
Filled Volume 100 ml
Storage Condition ≤ 8 degree Celsius or colder
Preservative(s) <Product Specific>
Filtration Product has been filtered through an absolute 0.2 micron filter membrane
Appearance Straw to light amber color with no color suggestive of hemolysis appearance
D.O.M. <Lot Specific>
Expiration Date 5 Years from the date of manufacture (D.O.M.)
Manufacturer Monobind Inc. ISO: 13485 Certified Company
Blood Testing All Units in this serum pool were found negative or non-reactive to HIV 1/2 Antibodies, HBsAg, HCV Antibodies, and HIV-Ag, HBV, and HCV by NAT, also ALT value less than the 2X upper limit of normal range by FDA approved methods.
Caution Not for use in human or clinical diagnosis. For use in research, or manufacturing non-injectable products only.

Analyte	Results	Methodology
Total Protein	TBD mg/dL	SMAC
Triglyceride	TBD mg/dL	SMAC
Chloride	TBD mg/dL	SMAC
Sodium	TBD IU/L	SMAC
Microbial Colony Count	TBD cfu/mL	Pour Plate, 30-35 degree C, 48 hrs.
pH	TBD	pH Meter

Prepared by _____ (Biomatrix Dept.)

Date: _____

Verified by _____ (QC Dept.)

Date: _____

