

CERTIFICATE OF ANALYSIS

Product Drug-Free Urine Matrix - Sodium Azide pH treated - 250 mL
Item No. 004CM-250ML
Lot No. 004CM-XXA24
Filled Volume 250 ml
Storage Condition ≤ 8 degree Celsius or colder
Preservative(s) <Product Specific>
D.O.M. <Lot Specific>
Expiration Date 5 Years from the date of manufacture (D.O.M.)
Manufacturer Monobind Inc. ISO: 13485 Certified Company
Donor Testing All Units in this urine pool were found negative or non-reactive for HIV 1/2 Antibodies, HBsAg, and HCV Antibodies by FDA approved methods.
Caution Not for use in human or clinical diagnosis. For use in research, or manufacturing non-injectable products only.
Urine GCMS/LCMS Testing No drugs detected

Test	Results	Methodology
Appearance	Clear Yellow	Visual
Uric Acid	TBD mg/dL	SMAC
pH	TBD	pH meter
Precipitation	No precipitation observed	Visual Observation
Microbial Colony Count	TBD cfu/mL	Pour Plate, 30-35 degree C, 48 hrs.

Prepared by _____ (Biomatrix Dept.)
Verified by _____ (QC Dept.)

Date: _____
Date: _____

