

100 North Pointe Drive Lake Forest, California 92630 USA Phone: 949-951-2665 Ex208 Fax: 949-951-3539

Certificate of Analysis

Product: Human Serum Specialty – Single Charcoal Stripped and Delipidated

Lot#: 198CM-01A09
Bulk: 1000 Liters
Filled in: 3.5 Liters/Bottle
pH: 7.45 (pH meter)

Preservative: None

Filtration: Filtered through a 0.22micron absolute membrane filter prior to bottling

Storage: -20°C

D.O.M.: 01-12-2009

Expiration date: 01-2014

Bleed testing: All units were tested at the donor level by FDA approved tests found

negative to antibodies HIV-I/II, HCV, HTLV-1, HBsAg and Syphilis.

Analyte	Specification and methodology	Result
Aldosterone	60 ng/dL by RIA	<3.0 ng/dL
α - Fetoprotein	150 ng/mL by IMx	1.4 ng/mL
CEA	25 ng/mL by IMx	<0.3 ng/mL
Cortisol	40 μg/dL by TDx	0.22 μg/dL
Estradiol	650 pg/mL by RIA	19.4 pg/mL
Ferritin	400 ng/mL by IMx	10 ng/mL
ß – hCG	150 mIU/mL by IMx	0.4 mIU/mL
hCG	150 mIU/mL by IMx	0.4 mIU/mL
Prolactin	40 ng/mL by IMx	0.0 ng/dL
PSA	10 ng/mL by IMx	<1.0 ng/dL
Insulin	120 μIU/mL by IMx	0.0 μIU/mL
PAP	50 ng/mL by IMx	0.4 ng/mL
Thyroglobulin	60 ng/mL by RIA	0.0 g/dL
Iron	140 µg/dL by Dimension	83 µg/dL
Ammonia	<300 µmol/L by SMAC	50 μmol/L
Total Protein	>3.0 g/dL by SMAC	4.5 g/dL

Caution: Not for use in products subject to license under section 351 of the public health safety acts. For use in research or manufacturing into non-injectable products only.