

CERTIFICATE OF ANALYSIS

Product Human Serum Albumin 5% Liquid Low Folate & Low B12 Reagent Grade
Master Lot No. 102RM-XXA24
Item No. HA05SLB-100ML
Filled Volume 100 mL
Storage Condition ≤ 8 degree Celsius or colder
Countries of Plasma Origin <Product Specific>
D.O.M. <Lot Specific>
Expiration Date 10 Years from the date of manufacture (D.O.M.)
Manufacturer Monobind Inc. ISO: 13485 : 2016 Certified Company
Bleed Testing All Units in this serum pool were found negative or non-reactive to HIV 1/2 Antibodies, HBsAg, HCV Antibodies, and HIV, HBV, and HCV by NAT, and Syphilis by FDA approved methods.
Fractionation Method Cohn-Oncley Cold Alcohol Procedure
Caution Not for use in human or clinical diagnosis. For use in research, or manufacturing non-injectable products only.
Note Results based on 5 % solution

Analyte	Results	Methodology
Purity	TBD %	Protein Electrophoresis
Major Analyte	TBD <Product Specific>	ELISA
Endotoxin	TBD EU/mL	SMAC
Moisture	TBD %	Thermogravimetric
Microbial	TBD cfu/mL	Pour Plate, 30-35 degree C, 48 hrs.

Prepared by _____ (Biomatrix Dept.)
Verified by _____ (QC Dept.)

Date: _____
Date: _____

