



MONOBIND, INC.
ISO 13485 Certified Company

Certificate of Analysis/Batch Release Certificate

Item#: 625MM

Monobind Lot #: 01G19

Product: 0 Calibrator Matrix (20.1ml) for Steroid (All) and Thyroid (T3 , T4)

Used for: (For Steroid (All) and Thyroid (T3 , T4)

Expiration Date: 2024-07-29

Date of Manufacture: 2019-07-29

Form: Clear yellow liquid filtered through an absolute 0.2µm membrane filter prior to bottling.

Storage: 2-8°C

Preservative: 0.02% Methylisothiazolone and 0.02% bromonitrodioxane

pH: 7.2 to 7.4

The results from tests for individual, proprietary components meet the specification determined for each of those components

Caution: Not for use in human or clinical diagnosis. For use in research or in manufacturing of non-injectable products only.

Authorized by: Janet Varona, Quality Control

Date: 2019-07-29