



AMH - Anti-Müllerian Hormone

AccuBind® ELISA & AccuLite® CLIA Microwells



Key Indicator

- Menopausal status
- Ovarian reserve



Other Uses

- Reproductive choices (freeze eggs/ try earlier)
- IVF prediction and drug dosing
- PCOS (high levels of AMH)
- Ovarian cancer treatment (GCT type)
- Ambiguous infant genitalia investigation



Monobind AMH Kit



INCUBATION TIME: 1 hour, no shaking
 ASSAY TIME: **1 hour, 20 minutes** (ELISA)
 SAMPLE SIZE: 50 µl
 SENSITIVITY: 0.04 ng/ml
 ELISA CORRELATION TO BECKMAN ACCESS:
R² 0.953 based on 142 patient samples

Expected Values in ng/ml

AGE RANGE	N	MEAN	MEDIAN
20-29	24	5.95	5.94
30-39	80	2.83	2.47
40-49	38	1.33	0.59

Monobind AMH AccuBind® ELISA
 Levels peak in puberty (1-15 ng/ml) and decline gradually until
 menopause (<1 ng/ml)

www.monobind.com

Plasma or Serum
 New Product Release: 4Q19





AccuBind® ELISA & AccuLite® CLIA Microwells

New Products

Bone Metabolism

- CALCITONIN
- PTH INTACT
- PTH WHOLE & INTACT VAST

Endocrine

- ACTH
- ALDOSTERONE

Steroids

- E1

Future Products

Fertility

- INHIBIN A
- INHIBIN B

Infectious

- TORCH (1Q20)

Steroids

- ANDROSTENEDIONE (4Q19)

www.monobind.com

TRY
Monobind Kits
TODAY



Tel +1.949.951.2665
Fax +1.949.951.3539
sales@monobind.com