

# Human Plasma Products

## Decision Guide (Item Number Selection)

### 1 Determine Your Base Material

|  | Source  | Recovered | Fresh Frozen | Off The Clot |
|--|---------|-----------|--------------|--------------|
| Protein Content (typical native)   | 5-5.3%  | 5.5-6%    | 5.5-6%       | 7%           |
| <b>2 Select The Product That Matches Your Needs:</b>                                       |         |           |              |              |
| <b>Gendered</b>  | Yes/No  | Yes/No    | Yes/No       | Yes/No       |
| <b>Unprocessed plasma (Pooled or Single Donor)</b>   | 173RM   | 177RM     | 174RM        | 245RM        |
| <b>Defibrinated Protease free</b>  | 174CM   | 018CM     | 170CM        | 245CM        |
| Protease free  | 200CM   | 204CM     |              |              |
| <b>Single Charcoal Stripped</b>  | 099CM   | 098CM     | 101CM        |              |
| Protease free  |         | 208CM     |              |              |
| <b>2X Charcoal Stripped</b>  | 103CM   | 100CM     | 104CM        |              |
| Protease free  |         | 209CM     |              |              |
| <b>Charcoal Stripped, Delipidated</b>  | 176CM   | 180CM     | 182CM        |              |
| Protease free  | 202CM   | 206CM     |              |              |
| <b>2X Charcoal Stripped, Delipidated</b>   | 177CM   | 190CM     | 191CM        |              |
| Protease free  | 203CM   | 207CM     |              |              |
| <b>Delipidated</b>   | 175CM   | 198CM     | 171CM        | 246CM        |
| Protease free  | 201CM   | 205CM     |              |              |
| <b>Resin Stripped</b>  | 173CM   | 097CM     | 096CM        |              |
| Protease free  |         |           |              |              |
| <b>Delipidated, Resin Stripped</b>   | 183CM   | 185CM     | 095CM        |              |
| Protease free  |         |           |              |              |
| <b>Speciality HCY Depleted, Low B12, Anti-Tg/TPO, 15% T-Uptake, Calibrator Matrix etc.</b> | Call Us | Call Us   | Call Us      | Call Us      |

### 3 Sodium Azide as a Preservative?

- Yes (We recommend 0.1% sodium azide)  
 No  
 Custom preservative (let us know)

### 4 Analyte Profile

- Total protein count  
 Specific analyte values (enzyme, calcium, ALP)  
 pH, electrolytes and other analytes ranges (provide targets)

### 5 Pack Size?

- 3.5 liter jugs (recommended to reduce freight/ ensure transport temp)  
 10 liter jugs  
 Other (1, 2, 3 liter or custom)

### 6 Any Other Special Requirements?

Let us know

- > Once an item is selected, please [Contact Us](#) with your required quantity so we may prepare a quote.
- > For small quantities, some items may be in stock and samples available (fee may apply).
- > For large quantities or custom jobs, Monobind will provide you a pre-shipment sample for functional testing.
- > Based on your feedback, the bulk is released or other actions determined.
- > All shipments contain product Certificate of Analysis. All other traceability records retained by Monobind for compliance purposes.